



PRECISION BALANCE BFT-1804

PRECISION BALANCE BFT-1804

PRECISION BALANCE



LCD/AC and DC power supply (AA*3)
 Tare function / Counting / Unit conversion (g/ct/oz...)
 Min weighing set
 Overload alarm / Level indicator
 Option: Dual display / Interface / Printer / Dust cover

SPECIFICATIONS

Model	BFT-1804
Capacity	110g/610g
Readability	0.001g/0.01g
Min Weighing	1d-4d
Repeatability	±2d
Linearity	±3d
Stable Time	≤3s
Pan Size	Φ110mm
Appearance Size	270*200*280mm
Draft Shield Size	170*165*180mm (Inner diameter)
Calibration	External calibration
Cal. Weight	200g
Packing Size	360*265*350mm / 560*380*730mm (4 units per box)
Gross Weight	4kg / 16kg (4 units per box)
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1mg-500mg

1



1g-200g

1



1kg-5kg

1



1



2



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com