



## PRECISION BALANCE BFT-1903

# PRECISION BALANCE BFT-1903

## PRECISION BALANCE



Aluminum shell / Large space wind proof cover  
LED (Green) / LCD  
Tare function / Counting  
Overload alarm / Level indicator  
Option: Under weighing / Printer

## SPECIFICATIONS

Model	BFT-1903
Capacity	310g
Readability	0.001g
Min Weighing	0.004g
Repeatability	±0.002g
Linearity	±0.003g
Stable Time	≤3s
Pan Size	Φ90mm
Appearance Size	365*210*355mm
Draft Shield Size	175*165*230mm (Inner diameter)
Interface	RS232 (Standard)
Calibration	External calibration
Cal. Weight	200g
Option	Printer/Under weighing
Packing Size	Inside:440*295*430mm Outside:490*335*485mm
Gross Weight	7.5kg
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1



1



1



2



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)