



## MOISTURE ANALYZER BFT-201

# MOISTURE ANALYZER BFT-201

BASIC MOISTURE ANALYZER



1. 5 Inch touch screen
2. Real time testing curve
3. Standard/Fast/Soft Mode
4. Automatic/Manual/Timed stop mode
5. Silent fast cooling fan
6. Adaptive input voltage
7. Test frequency selection
8. Clock function
9. Historical test report

## SPECIFICATIONS

Model	BFT-201
Capacity	110 g
Readability	0.01 g
Sensor	Load Cell
Pan Size	Φ100 mm (20pcs)
Display	5 Inch touch screen
Operate. Temp	5°C-35°C
Heating Source	Halogen lamp
Temp. Unit	PT-100-2
Temp. Range	40°C-160°C
Temp. Setting	1°C
Moisture. Range	10.00%-90.00%
Moisture Readability	0.10%
Testing Range	10.00%-90.00%
DRY% Readability	0.10%
Storage	20
Analysis Curve	Real time dynamic curve
Testing Height	30mm
LANGUAGE	EN/CN
Power supplier	220V 50Hz / 110V 60Hz
Dimension	340*215*170mm
Packing Size	415*310*255mm
G.Weight/N.Weight	4.0kg/3.5kg
Alt Name	Basic Moisture Analyzer

## ACCESSORIES FOR PURCHASE

No	Name	Size	Material	Mat
1	Aluminum sample tray 1	Outer circle Φ 90mm Inner circle Φ 83mm Plate height 4mm	Aluminum	50 pieces/box
2	Temperature kit			



1



2

## FEATURES

1. 5 Inch touch screen
2. Real time testing curve
3. Standard/Fast/Soft Mode
4. Automatic/Manual/Timed stop mode
5. Silent fast cooling fan
6. Adaptive input voltage
7. Test frequency selection
8. Clock function
9. Historical test report



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)