



PRECISION BALANCE BFT-2302

PRECISION BALANCE BFT-2302

PRECISION BALANCE



LCD display
 Stainless steel weighing pan (Both side can be used)
 USB charging port (No adapter)
 2 x AAA dry battery
 Units conversion: g / oz / ct / tl

SPECIFICATIONS

| | |
|-----------------|-------------------------------|
| Model | BFT-2302 |
| Capacity | 6000g |
| Readability | 0.5g |
| Repeatability | 1.0g |
| Linearity | 1.5g |
| Min Weighing | 2.0g |
| Stable Time | ≤3s |
| Pan Size | 157*157mm |
| Temperature | 5-35°C |
| Appearance Size | 160*225mm |
| Net weight | 0.6kg per unit / 15kg per box |
| Gross weight | 0.8kg per unit / 17kg per box |
| Alt Name | Precision Balance |

ACCESSORIES FOR PURCHASE

| No | Name | Description | Information |
|----|------------|---|--|
| 1 | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g |
| 2 | Weight set | Lock type series weights (5kg-20kg) | M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g |



1



1



1



1



2



2


APPLICATIONS



Baked Cake



Healthy Diet



Fruit Selection



Drug Weighing



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com