



## PRECISION BALANCE BFT-2303

# PRECISION BALANCE BFT-2303

PRECISION BALANCE



LCD display  
 Stainless steel weighing pan (Both side can be used)  
 USB charging port (No adapter)  
 2 x AAA dry battery  
 Units conversion: g / oz / ct / tl

## SPECIFICATIONS

Model	BFT-2303
Capacity	10000g
Readability	1g
Repeatability	2g
Linearity	3g
Min Weighing	4g
Stable Time	≤3s
Pan Size	157*157mm
Temperature	5-35°C
Appearance Size	160*225mm
Net weight	0.6kg per unit / 15kg per box
Gross weight	0.8kg per unit / 17kg per box
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1



1



1



2



2


## APPLICATIONS



Baked Cake



Healthy Diet



Fruit Selection



Drug Weighing



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)