



DENSITY BALANCE BFT-2402

DENSITY BALANCE BFT-2402

DENSITY BALANCE



Color plastic shell
LCD / Tare function
Overload alarm / Level indicator
Option: Interface (RS232 / RS485) / Printer
Solid / Liquid
Density frame (Solid / Liquid)
Direct display sample density
Density reference set
Density Components

SPECIFICATIONS

Model	BFT-2402
Type	Liquid
Capacity	310g
WR	0.01g
DR	0.001g/cm ³
Test Capacity	0.001-99.999g/cm ³
Pan Size-01	Φ50mm
Pan Size-02	90*90*25mm
Tank Size	50ml Beaker
Op-temp	5-35°C
Key name	Power,F,ENTER,CAL, 0
Cal. weight	200g
Material	Plastic shell
Outlook size	270*195*230mm
Packing Size	345*325*195mm
Weight	4.5kg
Density Kit	Liquid kit
Battery	Yes
Alt Name	Density Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g

2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g
---	------------	-------------------------------------	--



1



1



1



1



2



2

FEATURES

Color plastic shell

LCD / Tare function

Overload alarm / Level indicator

Option: Interface (RS232 / RS485) / Printer



Solid / Liquid

Density frame (Solid / Liquid)

Direct display sample density

Density reference set

Density Components



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com