



## DENSITY BALANCE BFT-2502

# DENSITY BALANCE BFT-2502

## DENSITY BALANCE



Color plastic shell / Aluminum shell  
LCD / Tare function  
Overload alarm / Level indicator  
Interface: RS232 (Standard)  
Option: RS485 / Printer / Wind shield  
Solid / Liquid / Dust / Absorbent  
Density frame (Solid / Liquid / Powder)  
Direct display sample density  
Density reference set

## SPECIFICATIONS

Model	BFT-2502
Type	Solid
Capacity	310g
WR	0.005g
DR	0.001g/cm <sup>3</sup>
Test Capacity	0.001-99.999g/cm <sup>3</sup>
Pan Size-01	60*152mm
Pan Size-02	70*110mm
Tank Size	139*100*78mm Water tank
Op-temp	5-35°C
Key name	Power,F,ENTER,CAL,0
Cal-weight	200g
Material	Aluminium shell
Outlook size	365*215*250mm
Packing Size	410*270*165mm
Weight	4.5kg
Density Kit	Solid kit
Glass Draft	Option
Battery	No
Alt Name	Density Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g

2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g
---	------------	-------------------------------------	--



1mg-500mg  
1



1g-200g  
1



1kg-5kg  
1



1



2



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)