



## DENSITY BALANCE BFT-2503

# DENSITY BALANCE BFT-2503

## DENSITY BALANCE



Color plastic shell / Aluminum shell  
LCD / Tare function  
Overload alarm / Level indicator  
Interface: RS232 (Standard)  
Option: RS485 / Printer / Wind shield  
Solid / Liquid / Dust / Absorbent  
Density frame (Solid / Liquid / Powder)  
Direct display sample density  
Density reference set

## SPECIFICATIONS

|               |                               |
|---------------|-------------------------------|
| Model         | BFT-2503                      |
| Type          | Solid                         |
| Capacity      | 610g                          |
| WR            | 0.01g                         |
| DR            | 0.001g/cm <sup>3</sup>        |
| Test Capacity | 0.001-99.999g/cm <sup>3</sup> |
| Pan Size-01   | 60*152mm                      |
| Pan Size-02   | 70*110mm                      |
| Tank Size     | 139*100*78mm Water tank       |
| Op-temp       | 5-35°C                        |
| Key name      | Power,F,ENTER,CAL, 0          |
| Cal-weight    | 200g                          |
| Material      | Aluminium shell               |
| Outlook size  | 365*215*250mm                 |
| Packing Size  | 410*270*165mm                 |
| Weight        | 4.5kg                         |
| Density Kit   | Solid kit                     |
| Glass Draft   | Option                        |
| Battery       | No                            |
| Alt Name      | Density Balance               |

## ACCESSORIES FOR PURCHASE

| No | Name       | Description                                             | Information                                                                                                                                          |
|----|------------|---------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 level!Suitable for 0.1g-0.01g<br>F2 level!Suitable for 0.01g-0.001g<br>F1 level!Suitable for 0.001g-0.0001g<br>E2 level!Suitable for over 0.0001g |

|   |            |                                     |                                                                                                             |
|---|------------|-------------------------------------|-------------------------------------------------------------------------------------------------------------|
| 2 | Weight set | Lock type series weights (5kg-20kg) | M1 levelSuitable for 0.1g-0.01g<br>F2 levelSuitable for 0.01g-0.001g<br>F1 levelSuitable for 0.001g-0.0001g |
|---|------------|-------------------------------------|-------------------------------------------------------------------------------------------------------------|



1



1




1



1



2



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)