



DENSITY BALANCE BFT-2509

DENSITY BALANCE BFT-2509

DENSITY BALANCE



Color plastic shell / Aluminum shell
LCD / Tare function
Overload alarm / Level indicator
Interface: RS232 (Standard)
Option: RS485 / Printer / Wind shield
Solid / Liquid / Dust / Absorbent
Density frame (Solid / Liquid / Powder)
Direct display sample density
Density reference set

SPECIFICATIONS

Model	BFT-2509
Type	Absorbent
Capacity	610g
WR	0.01g
DR	0.001g/cm ³
Test Capacity	0.001-99.999g/cm ³
Pan Size-01	60*152mm
Pan Size-02	70*110mm
Tank Size	139*100*78mm Water tank
Op-temp	5-35°C
Key name	Power,F,ENTER,CAL, 0
Cal-weight	200g
Material	Aluminium shell
Outlook size	365*215*250mm
Packing Size	410*270*165mm
Weight	4.5kg
Density Kit	Solid kit
Glass Draft	Option
Battery	No
Alt Name	Density Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level!Suitable for 0.1g-0.01g F2 level!Suitable for 0.01g-0.001g F1 level!Suitable for 0.001g-0.0001g E2 level!Suitable for over 0.0001g

2	Weight set	Lock type series weights (5kg-20kg)	M1 levelSuitable for 0.1g-0.01g F2 levelSuitable for 0.01g-0.001g F1 levelSuitable for 0.001g-0.0001g
---	------------	-------------------------------------	---



1



1




1



1



2



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com