



DENSITY BALANCE BFT-2601

DENSITY BALANCE BFT-2601

DENSITY BALANCE



Internal auto calibration / Electromagnetic force sensor
One key / Time set internal calibration
Aluminum shell / Large space wind proof cover
Tare function / Counting / Percentage
Under weighing / Unit conversion (g/mg/ct/oz...)
Sensitivity / Speed set / Overload alarm / Fault alarm / Level indicator
Density Components

SPECIFICATIONS

Model	BFT-2601
Type	solid
Capacity	220g
WR	0.0001g
DR	0.0001g/cm ³
Test Capacity	0.0001-99.9999g/cm ³
Pan Size-01	60*152mm
Pan Size-02	70*110mm
Tank Size	139*100*78mm WaterTank
Op-temp	15-25°C
Display	Touch panel
Cal. weight	Automatic internal calibration (default off) with 200g weight
Material	Aluminium shell
Outlook size	365*210*355 mm
Packing Size	465*315*660 mm
Weight	11 kg
Density Kit	Solid kit
Battery	No
Alt Name	Density Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 levelSuitable for 0.1g-0.01g F2 levelSuitable for 0.01g-0.001g F1 levelSuitable for 0.001g-0.0001g E2 levelSuitable for over 0.0001g

2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g
---	------------	-------------------------------------	--



1



1




1



1



2



2

FEATURES

Internal auto calibration / Electromagnetic force sensor

One key / Time set internal calibration



Aluminum shell / Large space wind proof cover

Tare function / Counting / Percentage

Under weighing / Unit conversion (g/mg/ct/oz...)

Sensitivity / Speed set / Overload alarm / Fault alarm / Level indicator

Density Components



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com