



DENSITY BALANCE BFT-2701

DENSITY BALANCE BFT-2701

HYDROSTATIC BALANCE



For highway and textile industry
 Metal base
 Large size stainless steel plates
 LED (Green)
 Tare function / Counting
 Overload alarm
 Option: Interface / Printer

SPECIFICATIONS

Model	BFT-2701
Capacity	5100g
Readability	0.1g
Min Weighing	0.4g
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	160*160mm
Temperature	5-35°C
Packing Size	485*340*350mm
Gross Weight	6kg
Alt Name	Hydrostatic Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1




1



1



2



2

APPLICATIONS



Road Construction Engineering



Textile industry



Colleges and University



Laboratory



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com