



TEXTILE BALANCE BFT-2801

TEXTILE BALANCE BFT-2801

RAW SILK SIZE TESTER



Aluminum shell / Large space wind proof cover
LCD (White back light with black font)
Tare function
Language: Chinese / English
Overload alarm / Level indicator
Standard: Interface RS232 / Printer
Date / Time / Number of samples set
Automatically store 200 data / Delete and continuation
Inquiry / Print: Testing date, time, checker, number of machinery, the distribution of raw silk size, hank number, sum, Average fineness, Denier deviation, maximum deviation.

SPECIFICATIONS

Model	BFT-2801
Capacity	110g (99.5D)
Readability	0.001g (0.5D)
Min Weighing	0.004g
Repeatability	±0.002g
Linearity	±0.003g
Pan Size	Φ90mm
Temperature	5-25°C
Interface	RS232
Draft Shield Size	175*165*230mm
Packing Size	490*335*485mm
Cal.Weight	100g
Gross Weight	9kg
Alt Name	Raw Silk Size Tester



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com