

TEXTILE BALANCE BFT-2901

COUNT BALANCE



Aluminum shell / Large space wind proof cover
LCD (White back light with black font)
Tare function
Available in English and Chinese
Overload alarm / Level indicator
Interface RS232 / Printer
Date / Time / Number of samples / Yarn blend codeset
Automatically store 100 data / Delete and continuation
Unit conversion: OZ / Tex / Nm / Ne / D / g / m² / oz / yd²
Statistics / Print: Each weighing value, Maximum value (Max), Minimum value (Min), Average (AV), Uneven rate (CV)

SPECIFICATIONS

Model	BFT-2901
Capacity	210g
Readability	0.001g
Min Weighing	0.004g
Repeatability	±0.002g
Linearity	±0.003g
Pan Size	Φ90mm
Temperature	5-25°C
Interface	RS232
Draft Shield Size	175*165*230mm
Packing Size	490*335*485mm
Cal.Weight	200g
Gross Weight	9kg
Alt Name	Count Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com