



TEXTILE BALANCE BFT-3304

TEXTILE BALANCE BFT-3304

HIGH CAPACITY BALANCES



- High precision electromagnetic force sensor
- Dot Matrix graphical display
- External calibration
- GLP/ISO record of last stored calibration data
- Multi Languages
- RS232 serial interface
- Formulation (Recipe)
- Check weighing
- Percentage weighing
- Accumulation function
- Animal weighing
- Max load determination (Peak hold)

SPECIFICATIONS

Model	BFT-3304
Capacity	12100g
Readability	0.01g
Pan Size	195*175mm
Repeatability	0.01g
Linearity	±0.03g
Stable Time	≤4S
Alt Name	High Capacity Balances



FEATURES

- High precision electromagnetic force sensor
- Dot Matrix graphical display
- External calibration
- GLP/ISO record of last stored calibration data
- Multi Languages
- RS232 serial interface
- Formulation (Recipe)
- Check weighing
- Percentage weighing
- Accumulation function
- Animal weighing
- Max load determination (Peak hold)
- Plastic protection cover
- Under weighing

APPLICATIONS



Food industry



Gold Jewellery



Scientific Research Experiment



Cloth Textile



Medicine



Parts Inspection



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com