



DENSITY BALANCE BFT1BB2

DENSITY BALANCE BFT1BB2

HYDROSTATIC BALANCE



For highway and textile industry
Metal base
Large size stainless steel plates
LED (Green)
Tare function / Counting
Overload alarm
Option: Interface / Printer

SPECIFICATIONS

Model	BFT1BB2
Capacity	6100g
Readability	0.1g
Min Weighing	0.4g
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	160*160mm
Temperature	5-35°C
Packing Size	485*340*350mm
Gross Weight	6kg
Alt Name	Hydrostatic Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1

APPLICATIONS



Road Construction Engineering



Textile industry



Colleges and University



Laboratory



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com