



TEXTILE BALANCE BFT1BC1

TEXTILE BALANCE BFT1BC1

RAW SILK SIZE TESTER



Aluminum shell / Large space wind proof cover

LCD (White back light with black font)

Tare function

Language: Chinese / English

Overload alarm / Level indicator

Standard: Interface RS232 / Printer

Date / Time / Number of samples set

Automatically store 200 data / Delete and continuation

Inquiry / Print: Testing date, time, checker, number of machinery, the distribution of raw silk size, hank number, sum, Average fineness, Denier deviation, maximum deviation.

SPECIFICATIONS

Model	BFT1BC1
Capacity	110g (99.5D)
Readability	0.001g (0.5D)
Min Weighing	0.004g
Repeatability	±0.002g
Linearity	±0.003g
Pan Size	Φ90mm
Temperature	5-25℃
Interface	RS232
Draft Shield Size	175*165*230mm
Packing Size	490*335*485mm
Cal.Weight	100g
Gross Weight	9kg
Alt Name	Raw Silk Size Tester



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com