



## TEXTILE BALANCE BFT1BD1

# TEXTILE BALANCE BFT1BD1

## COUNT BALANCE



Aluminum shell / Large space wind proof cover  
LCD (White back light with black font)  
Tare function  
Available in English and Chinese  
Overload alarm / Level indicator  
Interface RS232 / Printer  
Date / Time / Number of samples / Yarn blend codeset  
Automatically store 100 data / Delete and continuation  
Unit conversion: OZ / Tex / Nm / Ne / D / g / m<sup>2</sup> / oz / yd<sup>2</sup>  
Statistics / Print: Each weighing value, Maximum value (Max), Minimum value (Min), Average (AV), Uneven rate (CV)

## SPECIFICATIONS

| Model             | BFT1BD1       |
|-------------------|---------------|
| Capacity          | 210g          |
| Readability       | 0.001g        |
| Min Weighing      | 0.004g        |
| Repeatability     | ±0.002g       |
| Linearity         | ±0.003g       |
| Pan Size          | Φ90mm         |
| Temperature       | 5-25°C        |
| Interface         | RS232         |
| Draft Shield Size | 175*165*230mm |
| Packing Size      | 490*335*485mm |
| Cal.Weight        | 200g          |
| Gross Weight      | 9kg           |
| Alt Name          | Count Balance |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)