



TEXTILE BALANCE BFT1BE1

TEXTILE BALANCE BFT1BE1

DANIEL BALANCE



LED (Green)
Tare function / Counting
Overload alarm
Option: Under weighing / Interface / Printer
Number of samples set
Display: G (g) / Dan Neal value (D) (size value)

SPECIFICATIONS

Model	BFT1BE1
Capacity	110g (99.5D)
Readability	0.001g (0.5D)
Min Weighing	0.004g
Pan Size	Φ130mm
Temperature	5-25°C
Cal. Weight	100g
Packing Size	310*230*155mm
Gross Weight	3kg
Alt Name	Daniel Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com