



TEXTILE BALANCE BFT1BF1

TEXTILE BALANCE BFT1BF1

DANIEL BALANCE



LCD (White back light with black font)
AC and DC power supply (AA*3)
Tare function / Counting / Unit conversion (g, OZ, Tex, Nm, Ne, D, g/m², oz/yd²)
Overload alarm / Level indicator
Option: Dual display / Interface / Printer

SPECIFICATIONS

Model	BFT1BF1
Capacity	310g
Readability	0.01g
Min Weighing	0.04g
Pan Size	Φ130mm
Temperature	5-25°C
Cal. Weight	200g
Packing Size	310*230*155mm
Gross Weight	2.5kg
Alt Name	Daniel Balance

ACCESSORIES FOR PURCHASE

No	Name	Size
1	Plastic windshield	18*17*9.5cm
2	Collapsible plastic windshield	18*18.5*9cm
3	Dust cover	19*27*6.5cm



1



2



3



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com