



## TEXTILE BALANCE BFT1BH3

# TEXTILE BALANCE BFT1BH3

## HIGH CAPACITY BALANCES



High precision electromagnetic force sensor  
Dot Matrix graphical display  
External calibration  
GLP/ISO record of last stored calibration data  
Multi Languages  
RS232 serial interface  
Formulation (Recipe)  
Check weighing  
Percentage weighing  
Accumulation function  
Animal weighing  
Max load determination (Peak hold)

## SPECIFICATIONS

Model	BFT1BH3
Capacity	10100g
Readability	0.01g
Pan Size	195*175mm
Repeatability	0.01g
Linearity	±0.03g
Stable Time	≤4S
Alt Name	High Capacity Balances



## FEATURES

High precision electromagnetic force sensor  
Dot Matrix graphical display  
External calibration  
GLP/ISO record of last stored calibration data  
Multi Languages  
RS232 serial interface  
Formulation (Recipe)  
Check weighing  
Percentage weighing  
Accumulation function  
Animal weighing  
Max load determination (Peak hold)  
Plastic protection cover  
Under weighing

## APPLICATIONS



Food industry



Gold Jewellery



Scientific Research Experiment



Cloth Textile



Medicine



Parts Inspection



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)