



TEXTILE BALANCE BFT1BH4

TEXTILE BALANCE BFT1BH4

HIGH CAPACITY BALANCES



High precision electromagnetic force sensor
Dot Matrix graphical display
External calibration
GLP/ISO record of last stored calibration data
Multi Languages
RS232 serial interface
Formulation (Recipe)
Check weighing
Percentage weighing
Accumulation function
Animal weighing
Max load determination (Peak hold)

SPECIFICATIONS

Model	BFT1BH4
Capacity	12100g
Readability	0.01g
Pan Size	195*175mm
Repeatability	0.01g
Linearity	±0.03g
Stable Time	≤4S
Alt Name	High Capacity Balances



FEATURES

High precision electromagnetic force sensor
Dot Matrix graphical display
External calibration
GLP/ISO record of last stored calibration data
Multi Languages
RS232 serial interface
Formulation (Recipe)
Check weighing
Percentage weighing
Accumulation function
Animal weighing
Max load determination (Peak hold)
Plastic protection cover
Under weighing

APPLICATIONS



Food industry



Gold Jewellery



Scientific Research Experiment



Cloth Textile



Medicine



Parts Inspection



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com