



ACCESSORIES BFT1BJ1

ACCESSORIES BFT1BJ1

STATIC ELIMINATOR



SPECIFICATIONS

Model	BFT1BJ1
Distance with ion source and sample	5-40 cm
Discharge time	9 sec/5 cm, 13 sec/10 cm, 100 sec/40 cm (When fan is on)
Ozone concentration	0-0.05 ppm (Distance of ion source 2 cm)
Max. air volume	0.06 cm ³ /second
Environment condition	0-50°C, 20-80% Air Humidity (Non-condensing)
AC Adapter	AC 100-200V, 50-60Hz
Size	11*10.5*4.5 cm
Weight	310 g
Installation place	Indoors only
Alt Name	Static eliminator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com