



WEIGHING SCALE BFT1BL6 (BBAL-415)

WEIGHING SCALE BFT1BL6

INDUSTRIAL SCALE



Large range and high precision
Spray steel scale frame
Large size stainless steel table
LED (Red)
AC and DC power supply
Tare function / Counting / Unit conversion (kg/lb)
Overload alarm / Undervoltage alarm / Level indicator
Option: Interface / Printer

SPECIFICATIONS

Model	BFT1BL6
Old Model	BBAL-415
Capacity	510kg
Readability	50g
Repeatability	±100g
Linearity	±150g
Pan Size	600x800mm
Stable Time	≤2S
Alt Name	Industrial scale



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com