



## WEIGHING SCALE BFT1BL9 (BBAL-419)

# WEIGHING SCALE BFT1BL9

## INDUSTRIAL SCALE



Large range and high precision  
Spray steel scale frame  
Large size stainless steel table  
LED (Red)  
AC and DC power supply  
Tare function / Counting / Unit conversion (kg/lb)  
Overload alarm / Undervoltage alarm / Level indicator  
Option: Interface / Printer

## SPECIFICATIONS

|               |                  |
|---------------|------------------|
| Model         | BFT1BL9          |
| Old Model     | BBAL-419         |
| Capacity      | 160kg            |
| Readability   | 1g               |
| Repeatability | ±2g              |
| Linearity     | ±3g              |
| Pan Size      | 400x500mm        |
| Stable Time   | ≤2S              |
| Alt Name      | Industrial scale |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)