



WEIGHING SCALE BFT1BM3 (BBAL-408)

WEIGHING SCALE BFT1BM3

INDUSTRIAL SCALE



Large range and high precision
Fashionable appearance
Large size stainless steel pan
LCD (White back light with black font)
AC and DC power supply
Tare function / Counting /
Unit conversion (g/ct/oz), Option lb
Overload alarm / Level indicator
Option: Interface / Printer

SPECIFICATIONS

Model	BFT1BM3
Old Model	BBAL-408
Capacity	1g
Readability	21kg
Repeatability	±2g
Linearity	±3g
Pan Size	315x230mm
Stable Time	≤2S
Alt Name	Industrial scale



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com