



ANALYTICAL BALANCE BFT1K6

ANALYTICAL BALANCE BFT1K6

SEMI-MICRO ANALYTICAL BALANCE



Recipe database: 99 datas can be stored
 Automatic internal calibration
 Tare / counting / Percentage test / Animal weighing
 Baud rate is adjustable / Density function
 Unit conversion (g/mg/ct/oz...)
 Multi-Language: English / German / French / Italian / Spanish / Portuguese
 Keep weight values in GLP/ISO format, includes date, time, and serial number
 Plastic protection cover for keypad
 Overload alarm / Fault alarm / Level indicator
 Ioniser (option)

SPECIFICATIONS

Model	BFT1K6
Capacity	210g
Readability	0.01mg
Min weighing	100mg
Stable Time	≤6s
Display	LCD(white back light with black font)
Pan size	Φ80mm
Operation temp	5-35°C 20-80% Humidity
Repeat ability	0.05mg
Linearity	± 0.1mg
Draft shield size	152mmx160mmx225mm
Cal.weight	Internal calibration
Option	Printer
Dimension	345x215x345mm
G.N. weight	8.3kg/5.85kg
IONISER	Standard
Language	Italian, English, Portugues, Deutsch, France, Espanol
Packing Size	510x400x620mm
Power Supply	110-230V 50/60Hz Out DC 12V 1A
Alt Name	Semi-micro Analytical Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1

FEATURES

Recipe database: 99 datas can be stored

Automatic internal calibration

Tare / counting / Percentage test / Animal weighing

Baud rate is adjustable / Density function

Unit conversion (g/mg/ct/oz...)

Multi-Language: English / German / French / Italian / Spanish / Portuguese

Keep weight values in GLP/ISO format, includes date, time, and serial number

Plastic protection cover for keypad

Overload alarm / Fault alarm / Level indicator



Ioniser (option)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com