



## ANALYTICAL BALANCE BFT1N2

# ANALYTICAL BALANCE BFT1N2

INTERNAL CALIBRATION TOUCH SCREEN ANALYTICAL BALANCE



Internal auto calibration  
5 Inch touch screen  
One key / Time trigger internal calibration  
Electromagnetic force sensor  
Aluminum shell / Large space wind proof cover  
Tare function / Counting / Density / Percentage  
Unit conversion (g/mg/ct/oz...)  
Under weighing  
Sensitivity / Speed set  
Overload alarm / Fault alarm / Level indicator  
Option: Density kit (Solid / Liquid)

## SPECIFICATIONS

|                   |                                                      |
|-------------------|------------------------------------------------------|
| Model             | BFT1N2                                               |
| Capacity          | 110g                                                 |
| Readability       | 0.0001g                                              |
| Min Weighing      | 0.0004g                                              |
| Repeatability     | ±0.0002g                                             |
| Linearity         | ±0.0003g                                             |
| Stable Time       | ≤3s                                                  |
| Pan Size          | Φ80mm                                                |
| Appearance Size   | 365*210*355mm                                        |
| Draft Shield Size | 175*165*230mm (Inner diameter)                       |
| Communication     | RS232 (Standard) USB (Option)                        |
| Temperature       | 15-35°C                                              |
| Calibration       | Internal calibration                                 |
| Cal. Weight       | 100g                                                 |
| Option            | Printer                                              |
| Packing Size      | Inside:440*295*430mm Outside:490*335*485mm           |
| Gross Weight      | 11kg                                                 |
| Alt Name          | Internal Calibration Touch Screen Analytical Balance |

## ACCESSORIES FOR PURCHASE

| No | Name       | Description                                             | Information                      |
|----|------------|---------------------------------------------------------|----------------------------------|
| 1  | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 level Suitable for 0.1g-0.01g |



1



1





1



1

## FEATURES

Internal auto calibration  
 5 Inch touch screen  
 One key / Time trigger internal calibration  
 Electromagnetic force sensor  
 Aluminum shell / Large space wind proof cover  
 Tare function / Counting / Density / Percentage  
 Unit conversion (g/mg/ct/oz...)  
 Under weighing  
 Sensitivity / Speed set  
 Overload alarm / Fault alarm / Level indicator  
 Option: Density kit (Solid / Liquid)



## APPLICATIONS



Gold jewelry



High precision data matching



Scientific experiment



Medicine



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)