



ANALYTICAL BALANCE BFT1N4 (BBAL-111)

ANALYTICAL BALANCE BFT1N4

INTERNAL CALIBRATION TOUCH SCREEN ANALYTICAL BALANCE



Internal auto calibration
 5 Inch touch screen
 One key / Time trigger internal calibration
 Electromagnetic force sensor
 Aluminum shell / Large space wind proof cover
 Tare function / Counting / Density / Percentage
 Unit conversion (g/mg/ct/oz...)
 Under weighing
 Sensitivity / Speed set
 Overload alarm / Fault alarm / Level indicator
 Option: Density kit (Solid / Liquid)

SPECIFICATIONS

Model	BFT1N4
Old Model	BBAL-111
Capacity	200g
Readability	0.0001g
Min Weighing	0.0004g
Repeatability	±0.0002g
Linearity	±0.0003g
Stable Time	≤3s
Pan Size	Φ80mm
Appearance Size	365*210*355mm
Draft Shield Size	175*165*230mm (Inner diameter)
Communication	RS232 (Standard) USB (Option)
Temperature	15-35℃
Calibration	Internal calibration
Cal. Weight	200g
Option	Printer
Packing Size	Inside:440*295*430mm Outside:490*335*485mm
Gross Weight	11kg
Alt Name	Internal Calibration Touch Screen Analytical Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1





1



1

FEATURES

Internal auto calibration
5 Inch touch screen
One key / Time trigger internal calibration
Electromagnetic force sensor
Aluminum shell / Large space wind proof cover
Tare function / Counting / Density / Percentage
Unit conversion (g/mg/ct/oz...)
Under weighing
Sensitivity / Speed set
Overload alarm / Fault alarm / Level indicator
Option: Density kit (Solid / Liquid)



APPLICATIONS



Gold jewelry



High precision data matching



Scientific experiment



Medicine



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com