



ANALYTICAL BALANCE BFT103

ANALYTICAL BALANCE BFT103

INTERNAL CALIBRATION ANALYTICAL BALANCE



Internal auto calibration
 LCD (White back light with black font)
 One key / Time set internal calibration
 Electromagnetic force sensor
 Aluminum shell / Large space wind proof cover
 Tare function / Counting / Percentage
 Unit conversion (g/mg/ct/oz...)
 Under weighing
 Sensitivity / Speed set
 Overload alarm / Fault alarm / Level indicator

SPECIFICATIONS

Model	BFT103
Capacity	120g
Readability	0.0001g
Min Weighing	0.0004g
Repeatability	±0.0002g
Linearity	±0.0003g
Stable Time	≤3s
Pan Size	Φ80mm
Appearance Size	365*210*355mm
Draft Shield Size	175*165*230mm (Inner diameter)
Communication	RS232 (Standard) USB (Option)
Temperature	15-35°C
Calibration	Internal calibration
Cal. Weight	100g
Option	Printer
Packing Size	Inside:440*295*430mm Outside:490*335*485mm
Gross Weight	11kg
Alt Name	Internal Calibration Analytical Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com