



PRECISION BALANCE BFT1T1 (BBAL-101)

PRECISION BALANCE BFT1T1

PRECISION BALANCE



Aluminum shell / Large space wind proof cover
LED (Green) / LCD
Tare function / Counting
Overload alarm / Level indicator
Option: Under weighing / Printer

SPECIFICATIONS

Model	BFT1T1
Old Model	BBAL-101
Capacity	110g
Readability	0.001g
Min Weighing	0.004g
Repeatability	±0.002g
Linearity	±0.003g
Stable Time	≤3s
Pan Size	Φ90mm
Appearance Size	365*210*355mm
Draft Shield Size	175*165*230mm (Inner diameter)
Interface	RS232 (Standard)
Calibration	External calibration
Cal. Weight	100g
Option	Printer/Under weighing
Packing Size	Inside:440*295*430mm Outside:490*335*485mm
Gross Weight	7.5kg
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com