



PRECISION BALANCE BFT1U4 (BBAL-225)

PRECISION BALANCE BFT1U4

PRECISION BALANCE



Aluminum shell
Unit conversion (g/ct/oz...)
LED (Green) / LCD
Tare function / Counting
Overload alarm / Level indicator
Option: Printer

SPECIFICATIONS

Model	BFT1U4
Old Model	BBAL-225
Capacity	5000g
Readability	0.01g
Repeatability	±0.02g
Linearity	±0.03g
Pan Size	180*180mm
Stable Time	≤3s
Interface	RS232 (Standard)
Dimensions	345*205*105mm
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1

FEATURES



RS232

Level bubble

Power



Pan

Level bubble foot

LCD display screen

Aluminum shell
Unit conversion (g/ct/oz...)
LED (Green) / LCD
Tare function / Counting
Overload alarm / Level indicator
Option: Printer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com