



PRECISION BALANCE BFT1V1

PRECISION BALANCE BFT1V1

PRECISION BALANCE



LCD / AC and DC power supply (AA*3)
Tare function / Counting / Unit conversion (g/ct/oz...)
Min weighing set
Overload alarm / Level indicator
Option: Dual display / Interface / Printer / Windproof / Dust cover
Note: Capacity more than 1 kg or accuracy less than 0.01 g is not equipped with calibration weight.

SPECIFICATIONS

Model	BFT1V1
Capacity	110g
Readability	0.01g
Repeatability	±0.02g
Linearity	±0.03g
Pan Size	Φ130mm
Stable Time	≤2S
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1

No	Name	Size
1	Plastic windshield	18*17*9.5cm
2	Collapsible plastic windshield	18*18.5*9cm
3	Dust cover	19*27*6.5cm



1



2



3



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com