



PRECISION BALANCE BFT1W1

PRECISION BALANCE BFT1W1

PRECISION BALANCE



Fashional appearance
 Large size stainless steel pan (For XY-1BF)
 LCD (White back light with black font)
 AC and DC power supply (AA*3)
 Tare function / Counting / Unit conversion (g/ct/oz...)
 Min weighing set
 Overload alarm / Level indicator
 Option: Dual display / Interface / Printer / Dust cover
 Calibration weight

SPECIFICATIONS

Model	BFT1W1
Capacity	310g
Readability	0.1g
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	Φ130mm
Stable Time	≤2S
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Size
1	Plastic windshield	18*17*9.5cm
2	Collapsible plastic windshield	18*18.5*9cm
3	Dust cover	19*27*6.5cm



1



2



3

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1




1



1



2



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com