



PRECISION BALANCE BFT1X3

PRECISION BALANCE BFT1X3

PRECISION BALANCE



LCD display
Stainless steel weighing pan (Both side can be used)
USB charging port (No adapter)
2 x AAA dry battery
Units conversion: g / oz / ct / tl

SPECIFICATIONS

Model	BFT1X3
Capacity	10000g
Readability	1g
Repeatability	2g
Linearity	3g
Min Weighing	4g
Stable Time	≤3s
Pan Size	157*157mm
Temperature	5-35°C
Appearance Size	160*225mm
Net weight	0.6kg per unit / 15kg per box
Gross weight	0.8kg per unit / 17kg per box
Alt Name	Precision Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1

APPLICATIONS



Baked Cake



Healthy Diet



Fruit Selection



Drug Weighing



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com