



DENSITY BALANCE BFT1Z1 (BBAL-302)

DENSITY BALANCE BFT1Z1

DENSITY BALANCE



Color plastic shell / Aluminum shell
 LCD / Tare function
 Overload alarm / Level indicator
 Interface: RS232 (Standard)
 Option: RS485 / Printer / Wind shield
 Solid / Liquid / Dust / Absorbent
 Density frame (Solid / Liquid / Powder)
 Direct display sample density
 Density reference set

SPECIFICATIONS

Model	BFT1Z1
Old Model	BBAL-302
Type	Solid
Capacity	210g
WR	0.001g
DR	0.001g/cm ³
Test Capacity	0.001-99.999g/cm ³
Pan Size-01	60*152mm
Pan Size-02	70*110mm
Tank Size	139*100*78mm Water tank
Op-temp	5-35°C
Key name	Power,F,ENTER,CAL, 0
Cal-weight	200g
Material	Aluminium shell
Outlook size	365*215*250mm
Packing Size	410*270*165mm
Weight	4.5kg
Density Kit	Solid kit
Glass Draft	Option
Battery	No
Alt Name	Density Balance

ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g



1



1



1



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com