



CONDUCTIVITY METER BFU-2503

CONDUCTIVITY METER BFU-2503

PEN TYPE PH/CONDUCTIVITY METER



1. Mini design with size 180x50x20mm.
2. Auto calibration, auto power-off, low voltage alarm, electrode failure warning.
3. IP57 waterproof, replaceable electrode.

SPECIFICATIONS

Model	BFU-2503
Description	Conductivity Meter
Measuring Range	Conductivity (0 ~ 10)mS/cm; (0~99.9) μ S/cm, (100~999) μ S/cm, (1.00~10.00)mS/cm.
Resolution	0.1/1 μ S/cm; 0.01mS/cm
Accuracy	\pm 3%FS
Temp. Compensation	0 ~ 60°C
Calibration	1 Point
Power	Lithium batteries (CR2032) x 2
Humidity	0.86
Standard Accessories	Backup batteries
Alt Name	Pen Type pH/Conductivity Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com