



## PH METER BFU-2601

# PH METER BFU-2601

PH METER



1. Large LCD screen with controllable white backlight.
2. Built-in microprocessor chip has Intelligent functions such as auto calibration, ATC, data storage, function settings, auto diagnosis, auto power off, and undervoltage prompt, etc.
3. Auto recognition of 25 buffers with three kinds of options: Europe & USA, NIST and China. Support 1, 2, 3 point calibration.

## SPECIFICATIONS

Model	BFU-2601
pH Measuring range	(-2.00~19.99)pH
pH Resolution	0.1/0.01 pH
pH Accuracy	Meter:±0.01pH; Overall: ±0.02pH
pH Input impedance	≥1x10 <sup>12</sup> Ω
pH Stability	±0.01 pH/3h
mV Measuring range	±1999.9mV
mV Resolution	1mV
mV Accuracy	±0.1% FS
Temp. Measuring range	(0 ~ 100)°C
Temp. Resolution	0.1°C
Temp. Accuracy	5 ~ 60°C: ±0.4°C Others: ±0.8°C
Other Data storage	200 sets
Other Power	AA batteries (1.5V x2)
Other Ambient Temp.	5 ~ 35 °C
Other Ambient humidity	≤85%
Other IP grade	IP57 Dustproof and waterproof
Other Dimensions/N.W.	160x90x30mm/360g
Other Shipping size/G.W.	480x300x60mm/1KG
Standard Accessories	pH electrode, Temperature probe, pH buffer solution, Manual
Alt Name	pH Meter



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)