



## PH METER BFU-5202

# PH METER BFU-5202

PH METER



## SPECIFICATIONS

Model	BFU-5202
Parameters	pH/mV/°C
Range pH	-2 ~ 18.00pH
Range mV	±1999mV
Range Temp	(0.0 ~ 100)°C
Resolution Ph	0.01 pH
Resolution mV	1mV
Resolution Temp	0.1°C
Accuracy pH	±0.01pH±1bit pH
Accuracy mV	±0.1%FS
Accuracy Temp	±0.3°C±1bit
Temp. Compensation	(0.0 ~ 100)°C, auto/manual
Calibration	2 point
Power	220V/110V
Dimension/N.W.	300x200x90 mm / 1.5kg
Shipping size/G.W.	330x260x280 mm / 2kg
Standard Accessories	Electrode Holder
Alt Name	pH Meter



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)