



PH METER BFU-5204

PH METER BFU-5204

PH METER



SPECIFICATIONS

Model	BFU-5204
Parameters	Conductivity/TDS/°C
Range Cond	1.0.00 μ S/cm-100mS/cm
Range TDS	0-1000 mg/L
Range Temp	(0.0 ~ 60.0)°C
Accuracy Cond	\pm 0.5%FS
Accuracy TDS	\pm 1.0%FS
Accuracy Temp	\pm 0.3°C \pm 1bit
Temp. Compensation	(0.0 ~ 40.0)°C, auto/manual
Power	220V/110V
Dimension/N.W.	300x200x90 mm / 1.5kg
Shipping size/G.W.	
Standard Accessories	Electrode Holder
Alt Name	Conductivity meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com