



ANALYTICAL BALANCE BFU-5507

ANALYTICAL BALANCE BFU-5507

TOUCH SCREEN BALANCE



1. Large touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Standard anti-static module
4. 0.01mg maximum weighing range 100g
5. Electronic level adjustment indicator module
6. GLP/GMP

SPECIFICATIONS

Model	BFU-5507
Weighing Range	0-50g/0-100g
Calibration	External
Readability	0.01/0.1mg
Stabilization Time	0.01mg≤13S, 0.1mg≤4S
Pan Size	Φ80mm
Dimension/ N.W.	360*220*350mm/8KG
Shipping Size/ G.W.	520*380*500mm/10kg
Alt Name	Touch Screen Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com