



PRECISION BALANCE BFU-5608

PRECISION BALANCE BFU-5608

TOUCH SCREEN ELECTRONIC BALANCE



1. 7 inches touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing.
3. GLP/GMP

SPECIFICATIONS

| | |
|----------------------|--|
| Model | BFU-5608 |
| Weighing Range | 0~4200g |
| Readability | 10mg |
| Pan Size | 255*190mm |
| Dimension/N.W. | 280*200*60mm/ 1.8kg |
| Shipping Size/G.W. | 320*280*150mm/ 2.5kg |
| Standard Accessories | Calibration Weight (optional), Adapter, Manual, Dust cover |
| Alt Name | Touch Screen Electronic Balance |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com