



ANALYTICAL BALANCE BFU-5918

ANALYTICAL BALANCE BFU-5918

TOUCH SCREEN BALANCE



1. 4.3-inch color touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Electromagnetic balance sensor.
4. External/ Internal calibration.
5. Supports PC connection. Wireless bluetooth printing is optional.

SPECIFICATIONS

| | |
|--------------------|----------------------|
| Model | BFU-5918 |
| Maximum Capacity | 260g |
| Calibration | Internal |
| Readability | 1mg |
| Stabilization Time | ≤4S |
| Pan Size | Φ85mm |
| Dimension/N.W. | 360*220*350mm / 8kg |
| Shipping Size/G.W. | 520*380*500mm / 10kg |
| Alt Name | Touch Screen Balance |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com