



ANALYTICAL BALANCE BFU-6010

ANALYTICAL BALANCE BFU-6010

BALANCE



1. Black backlit 120*55 mm LCD screen
2. Electromagnetic balance sensor
3. External/ Internal calibration
4. 4-level shockproof
5. Counting function
6. Supports PC connection. Wireless bluetooth printing is optional.

SPECIFICATIONS

Model	BFU-6010
Weighing Range	0-210g
Calibration	Internal
Readability	0.1mg
Pan Size	Φ85mm
Dimension/N.W.	350*215*340mm/8KG
Shipping Size/G.W.	520*380*500mm/10kg
Standard Accessories	Calibration Weight (for external calibration), Adapter, Manual, Dust cover
Alt Name	0.1mg Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com