



## ANALYTICAL BALANCE BFU-6011

# ANALYTICAL BALANCE BFU-6011

BALANCE



1. Black backlit 120\*55 mm LCD screen
2. Electromagnetic balance sensor
3. External/ Internal calibration
4. 4-level shockproof
5. Counting function
6. Supports PC connection. Wireless bluetooth printing is optional.

## SPECIFICATIONS

Model	BFU-6011
Weighing Range	0-220g
Calibration	Internal
Readability	0.1mg
Pan Size	Φ85mm
Dimension/N.W.	350*215*340mm/8KG
Shipping Size/G.W.	520*380*500mm/10kg
Standard Accessories	Calibration Weight (for external calibration), Adapter, Manual, Dust cover
Alt Name	0.1mg Balance



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)