



## ANALYTICAL BALANCE BFU-6103

# ANALYTICAL BALANCE BFU-6103

BALANCE



1. Black backlit 120\*55 mm LCD screen
2. Electromagnetic balance sensor
3. External/ Internal calibration
4. 4-level shockproof
5. Counting function
6. Supports PC connection. Wireless bluetooth printing is optional

## SPECIFICATIONS

|                      |                                                                                    |
|----------------------|------------------------------------------------------------------------------------|
| Model                | BFU-6103                                                                           |
| Weighing Range       | 0-200g                                                                             |
| Calibration          | External                                                                           |
| Readability          | 1mg                                                                                |
| Pan Size             | Φ85mm                                                                              |
| Dimension/N.W.       | 350*215*340mm/8KG                                                                  |
| Shipping Size/G.W.   | 520*380*500mm/10kg                                                                 |
| Standard Accessories | Calibration Weight (for external calibration balance), Adapter, Manual, Dust cover |
| Alt Name             | Balance                                                                            |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)