



ANALYTICAL BALANCE BFU-6108

ANALYTICAL BALANCE BFU-6108

BALANCE



1. Black backlit 120*55 mm LCD screen
2. Electromagnetic balance sensor
3. External/ Internal calibration
4. 4-level shockproof
5. Counting function
6. Supports PC connection. Wireless bluetooth printing is optional

SPECIFICATIONS

Model	BFU-6108
Weighing Range	0-500g
Calibration	External
Readability	1mg
Pan Size	Φ85mm
Dimension/N.W.	350*215*340mm/8KG
Shipping Size/G.W.	520*380*500mm/10kg
Standard Accessories	Calibration Weight (for external calibration balance), Adapter, Manual, Dust cover
Alt Name	Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com