



ANALYTICAL BALANCE BFU-6109

ANALYTICAL BALANCE BFU-6109

BALANCE



1. Black backlit 120*55 mm LCD screen
2. Electromagnetic balance sensor
3. External/ Internal calibration
4. 4-level shockproof
5. Counting function
6. Supports PC connection. Wireless bluetooth printing is optional

SPECIFICATIONS

Model	BFU-6109
Weighing Range	0-100g
Calibration	Internal
Readability	1mg
Pan Size	Φ85mm
Dimension/N.W.	350*215*340mm/8KG
Shipping Size/G.W.	520*380*500mm/10kg
Standard Accessories	Calibration Weight (for external calibration balance), Adapter, Manual, Dust cover
Alt Name	Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com