



PH METER BFU1CA1

PH METER BFU1CA1

PH METER



SPECIFICATIONS

Model	BFU1CA1
Parameters	pH/mV
Range pH	0 ~ 14.00pH
Range mV	±1999mV
Range Temp	-
Resolution Ph	0.01 pH
Resolution mV	1mV
Resolution Temp	-
Accuracy pH	±0.01pH±1bit pH
Accuracy mV	±0.1%FS
Accuracy Temp	-
Temp. Compensation	(0.0~ 60.0)°C, manual
Calibration	1 or 2 point
Power	220V/110V
Dimension/N.W.	300x200x90 mm / 1.5kg
Shipping size/G.W.	330x260x280 mm / 2kg
Standard Accessories	Electrode Holder
Alt Name	pH Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com