



## PH METER BFU1CA4

# PH METER BFU1CA4

PH METER



## SPECIFICATIONS

|                      |                             |
|----------------------|-----------------------------|
| Model                | BFU1CA4                     |
| Parameters           | Conductivity/TDS/°C         |
| Range Cond           | 1.0.00μS/cm-100mS/cm        |
| Range TDS            | 0-1000 mg/L                 |
| Range Temp           | (0.0 ~ 60.0)°C              |
| Accuracy Cond        | ±0.5%FS                     |
| Accuracy TDS         | ±1.0%FS                     |
| Accuracy Temp        | ±0.3°C±1bit                 |
| Temp. Compensation   | (0.0 ~ 40.0)°C, auto/manual |
| Power                | 220V/110V                   |
| Dimension/N.W.       | 300x200x90 mm / 1.5kg       |
| Shipping size/G.W.   |                             |
| Standard Accessories | Electrode Holder            |
| Alt Name             | Conductivity meter          |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)