



## ANALYTICAL BALANCE BFU1CD11

# ANALYTICAL BALANCE BFU1CD11

## TOUCH SCREEN BALANCE



1. Large touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Standard anti-static module
4. 0.01mg maximum weighing range 100g
5. Electronic level adjustment indicator module
6. GLP/GMP

## SPECIFICATIONS

|                     |                      |
|---------------------|----------------------|
| Model               | BFU1CD11             |
| Weighing Range      | 200g                 |
| Calibration         | External             |
| Readability         | 0.1mg                |
| Stabilization Time  | ≤2S or 4S Switchable |
| Pan Size            | Φ80mm                |
| Dimension/ N.W.     | 360*220*350mm/8KG    |
| Shipping Size/ G.W. | 520*380*500mm/10kg   |
| Alt Name            | Touch Screen Balance |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)